

# Server Installation




Click [here](#) for system requirements.


---

## Required Prerequisites

1. **HTTP Activation must be turned on for the server.** This can be done in Server Manager's Roles and Features section under .Net Framework 3.5.
2. If installing on a PC make sure it does not go into power saving mode or sleep.

### Step OneStep TwoStep ThreeStep Four

	Click <a href="#">here</a> to download SureTrend and run it once completed.
	When the screen appears choose a setup language and click <b>Next</b> .
	Read and accept the license agreement and click <b>Next</b> to continue.

	<p>When you get to the <b>Multiuser - Distributed Network</b> checkbox check it. Just click <b>Next</b> to continue on with the installation.</p>
-----------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------


## Pointing clients to the server

Step OneStep TwoStep Three

Install SureTrend on any machine you wish to connect to the server **without checking the Multiuser - Distributed Network checkbox during installation.**

Start SureTrend on the client machine once installed and click on the **Admin** tab across the top. From the list of tabs on the left choose **Server Configuration Setting**.

Change the **Config Type** to Distributed and enter the Hostname or IP address of the server into the **Server** box. Click **Save** to restart SureTrend and when you've reopened it verify it is connected to the server by revisiting this area. *As shown below in Figure A.*



The image shows a configuration form for Hygiena, enclosed in a large blue circle. The form contains the following elements:

- Config Type:** A dropdown menu with "Distributed" selected.
- Server:** A text input field containing the placeholder text "IP Address or Hostname of Server".
- Port:** A text input field containing the value "83".
- Save:** A button located at the bottom of the form.

**Contact Hygiena Technical Support for further assistance.**

- Phone: 1-888-HYGIENA (1-888-494-4362, option 2)
- Email: [techsupport@hygiena.com](mailto:techsupport@hygiena.com)
- [Submit a Support Ticket](#)
- [Schedule a Microsoft Teams meeting with support](#)