


# Optimal Power Management Settings for your EnSURE™ Touch

## Optimal Display and Power Settings

1. From the **Home Screen**, tap **Settings**.
2. Tap **General**.
3. Tap **Display & Power Off**.
4. Adjust the **Display & Power Off** settings as listed below.
  - Set Display **Brightness** to 40% or less.
  - Set Display **Timeout** to 1-5 minutes.
  - Set Power Off **After Display Timeout** to **Off**.
  - Set Power Off **Daily at** to **Off**.




## Power Saving Recommendations

- Letting the screen sleep between tests or putting the screen to sleep by pressing the power button once will extend the battery charge.
- Turn the Wi-Fi off while the device is in use can help to increase the

battery life on your luminometer. You only need the Wi-Fi on when syncing, updating the software, or getting remote technical support.

- Go to Settings > Wi-Fi and tap the toggle to turn Wi-Fi off.
- Make certain that the EnSURE Touch is properly powered off after use by holding the power for about 1-2 seconds then selecting **Power Off** from the menu that appears.



- Enabling the 'After Display Timeout' and/or the 'Daily At' will result in your luminometer automating the shutdown process thus requiring a power reset sequence to reinitiate the operating system. If pressing the power button does not initiate the operating system, please try a 12 second hold on the power button.

**If you are having trouble with the battery life on your EnSURE Touch, please contact Hygiena Technical Support for assistance.**

- Phone: 1-888-HYGIENA (1-888-494-4362), option 2
- Email: [techsupport@hygiena.com](mailto:techsupport@hygiena.com)
- [Submit a Support Ticket](#)
- [Schedule a Microsoft Teams meeting](#)