




Server Installation


Click [here](#) for system requirements.

Required Prerequisites

1. **HTTP Activation must be turned on for the server.** This can be done in Server Manager's Roles and Features section under .Net Framework 3.5.
2. If installing on a PC make sure it does not go into power saving mode or sleep.

Step One Step Two Step Three Step Four

 The image shows the SureTrend logo, which consists of a blue circle containing a white square with a blue and yellow stylized 'S' and 'T' inside. Below the square, the word 'SureTrend' is written in white.	<p>Click here to download SureTrend and run it once completed.</p>
 The image shows a language selection screen with a list of languages: Chinese (Traditional), Czech, Dutch, English (United States), French (Standard), Italian, Portuguese (Standard), Russian, and Spanish. The 'English (United States)' option is highlighted with a blue bar.	<p>When the screen appears choose a setup language and click Next.</p>
 The image shows a license agreement screen with two radio buttons. The first option is 'I accept the terms of the license' and the second is 'I do not accept the terms of the license'. The first option is selected. There is a '< Back' button at the bottom.	<p>Read and accept the license agreement and click Next to continue.</p>



Pointing clients to the server

Step One Step Two Step Three

Install SureTrend on any machine you wish to connect to the server **without checking the Multiuser - Distributed Network checkbox during installation.**

Start SureTrend on the client machine once installed and click on the **Admin** tab across the top. From the list of tabs on the left choose **Server Configuration Setting.**

Change the **Config Type** to Distributed and enter the Hostname or IP address of the server into the **Server** box. Click **Save** to restart SureTrend and when you've reopened it verify it is connected to the server by revisiting this area. *As shown below in Figure A.*

The image shows a configuration form for Hygiena software, enclosed in a large blue circle. The form contains the following fields and controls:

- Config Type:** A dropdown menu with the value "Distributed" selected.
- Server:** A text input field containing the placeholder text "IP Address or Hostname of Server".
- Port:** A text input field containing the value "83".
- Save:** A button labeled "Save" located below the other fields.

Contact Hygiena Technical Support for further assistance.

- Phone: 1-888-HYGIENA (1-888-494-4362, option 2)
- Email: techsupport@hygiena.com
- [Submit a Support Ticket](#)
- [Schedule a Microsoft Teams meeting with support](#)