


Report Preview



1 - Search

Allows you to search for any text present in the report. This is useful for finding information based on specific keywords, phrases, etc.

2 - Open

Open any saved Preview Document Files (.prnx).

3 - Save

Save any preview to a Preview Document File (.prnx) on your computer.

4 - Print

Open your print configuration screen to print the report and modify print settings.

5/6 - Scale/Zoom

Adjust scaling of the report to fit perfectly on your page. Zoom in or out or select a specific zoom setting to view the report easier.

7 - Export

Export the report to a file. Choose from various document types like PDF, CSV, XLS or TXT; plus modify export settings.

8 - Send

Export and email the report. Choose from various document types like PDF, CSV, XLS, or TXT; plus modify export settings before sending.

9 - Watermark

Create a textual or picture watermark and apply it to the report.