

Download Center

IMPORTANT NOTICE ON UPDATING SURETREND

If you plan on updating your version of SureTrend to the latest, performing a backup first is essential to prevent lost or missing data. Follow the update instructions below and [contact us](#) if you need further assistance.

[Mandatory Update Instructions](#)

Installing this application requires administrative privileges. You may need your IT teams support if you are unable to install software on your machine.

This software is not intended for use with an EnSURE Touch system and is not necessary for use with SureTrend Cloud.

SureTrend is freeware. This means it is free software and you may install it freely on any compatible machine running Windows.

This software cannot be installed on Android phones or tablets, iPhones, on Macs, or any cloud-based services such as Drop Box, Google Drive, SharePoint, OneDrive, etc.

- **Operating System** = Windows XP and Windows Server 2008 R2 or newer
- **Memory** = 2 GB recommended
- **Disk Space** = 2 GB recommended
- **USB Requirements** = Must use supplied cable to sync

Choose an installation type below

- ☐

Standard Installation

Click [here](#) to view instructions on how to install SureTrend on one PC.

- ☐

Server Installation

Click [here](#) to view instructions on how to install SureTrend on a server and connect to it to share data.

- ☐

File Share Installation

Click [here](#) to view instructions on how to install SureTrend on one PC, move its database file to a network file share and connect to it to share its data.

- ☐

Updating SureTrend to the Latest Version

Click [here](#) to view instructions on how to update SureTrend and minimize the risk of losing your data.

Contact Hygiena Technical Support for further assistance.

- Phone: 1-888-HYGIENA (1-888-494-4362, option 2)
- Email: techsupport@hygiena.com
- [Submit a Support Ticket](#)
- [Schedule a Microsoft Teams meeting with support](#)

Download Details

Release date:	September 7, 2018
Last updated:	January 9, 2019
Current version:	4.0.3.4 - Build 11
Product type:	Software
File format:	.exe
File size:	168.33 MB

[Download Now](#)