

Results



The screenshot shows a mobile application interface for viewing results. At the top, there is a blue header with a back arrow and the word "Results". Below the header are four filter fields: "From" (text input), "To" (text input), "Plan" (dropdown menu), and "Location" (dropdown menu). Below the filters is a table with three columns: "Results", "%", and "Count 5". The table has three rows: "Pass" (green background, 60%, 3), "Caution" (yellow background, 0%, 0), and "Fail" (red background, 40%, 2). Each row has a right-pointing arrow next to the count.

Results	%	Count 5
✓ Pass	60%	3
⚠ Caution	0%	0
✗ Fail	40%	2

[Click to enlarge](#)

The Results screen will display all the results currently on your unit. You will be able to filter by date or choose to only show results from a specific plan or location.

The table area below the filters will display a total count of all Pass, Caution, and Fail results based on the filters selected above.